

## PROPRIO.springi

soft pad – plastic shell blank with sensomotoric POWERSPOTS



- infospots activate the musculature
- dimensionally stable bed ensures good heel hold
- flexible edges give pressure weals no chance
- certified material quality

**Express service: model 5130**

model	design	width	sizes
5130 long sole	long sole, blue microfibre cover	medium	19/20 – 35/36
5330 langsohlig	long sole, coloured textile cover- little bear	medium	19/20 – 35/36

processing – advantage	marketing-service
Position and height of the sensomotoric pulsers are changeable easily by spot warming. With PROPRIO.springi® you as a kid's specialist finally have the opportunity to treat the infantile talipes valgus cost- effectively according to the latest scientific findings.	<ul style="list-style-type: none"> <li>• flyers, deco column shop, window decoration</li> </ul>

active principle
<ul style="list-style-type: none"> <li>• support under the heel bone balcony for more intensive activation of the shin musculature</li> <li>• stabilisation of the foot contour by more intensive activation of the der calf bone musculature</li> <li>• flat floor contact of the toe tips</li> </ul>