

INJECTION MOULDED PLASTIC CORES

injection moulded plastic cores, stable, with slight heel hollow, lateral guide and strongly marked longitudinal arch



- in double sizes
- very good supporting core
- flexible heel
- free metatarsal bones
- super-thin

Express service: models 1802-00, 1803-00 and 1804-00

model	design	width	sizes
1802-00	half sole	slim	35/36 - 43/44
1803-00	half sole	medium	35/36 - 47/48
1804-00	half sole	wide	37/38 - 47/48

processing-advantages
<ul style="list-style-type: none"> • thermoformable at 160°-180°C • easy to grind with 50' or 60' grain at 1500 revolutions/min • glue advice: PP-glue e.g. ORTEC from Renia

additional design
<ul style="list-style-type: none"> • loose pad • lining: paraskin, loose/tight in skin or black • lining: microcork, 2 mm tight